pro 3d ure medical

AUDIO

24 | 25

Inhalt

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Contents

2 Introduction
4 **OTOFLOW® approved!**
6 **MSI® technology**
10 **audioprint®** 3D printing materials
18 Lightcuring gels
19 **otosil®** silicone materials
24 **SysTherm®** product family
26 Traditional production | cleaning
28 Accessories
30 Devices (scanner, software, 3D printer, cleaning, post curing)

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Materials are our DNA! ... is the guiding principle of pro3dure medical. Just like in our name, the 3d is the central pivot of all our activities and has been since 2013. With the purchase of land for the new headquarters of pro3dure, we will build a competence center for medical 3D printing in the Dortmund High Technology Park in the next few years in order to take the next step for you in terms of highly innovative material concepts. Keep an eye on our construction diary, which will appear on our website soon. **(Fig. 1–2)**

This year's focus on audiology is the result of our collaboration with Penrhos Bio, a subsidiary of Unilever. Using the globally unique **MSI® technology** (Multi Species Inhibition), biomimetic polymers for additive and traditional procedures that reduce the adhesion and formation of biofilms are being demonstrated for the first time. Accordingly, our materials are not only constantly improving in terms of mechanical properties, but are also becoming smarter in order to achieve better results. We have set ourselves the goal of learning from and with you in order to innovatively advance the manufacture of earmolds in all areas. Trusting cooperation and proximity to our customers and cooperation partners is therefore of outstanding importance to us in our daily work. Let's take the next steps together!

Abb. 1: Das zukünftige pro3dure Headquarter im Hochtechnologie-Park Dortmund.
Fig. 1: The future pro3dure headquarter in the Dortmund high-tech park.

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The flexible **OTOFLOW®** concept paves the way to a validated process.

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The demands for validated processes on the part of the approval authorities will increase for all market participants in the future. **pro3dure** has developed the **OTOFLOW® concept**, which can be used to compile a validated overall manufacturing process from verified process steps. This holistic approach ensures that not only the interactions between different materials, but also the individual process steps such as cleaning, post-curing and finishing are taken into account. The question **"Am I building the right earmold?"** can only be answered with a clear **"Yes"** in the future by taking an overall view.

MSI® technology

Inspired by Nature | Addicted to Function

Scotland's Life Sciences
Annual Awards
2024

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Bacteria have a fascinating and diverse social life. They exhibit coordinated group behavior controlled by **quorum-sensing systems** that detect the density of other bacteria in their environment. The term quorum dates back to the Roman Empire and referred to the lowest number of members required for a vote in the Senate. Similar to the Roman Senate, bacteria vote their behavior as soon as a certain number of other bacteria are in their vicinity. A key example of such group behavior is biofilm formation, in which communities of cells attach to a surface and form a polymer shell that protects the bacteria from external influences.

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Communication is everything and represents one of the most fundamental mechanisms of intra- and interspecies interactions. Communication can take place in many different ways. One form of human communication is the acoustic transmission of signals to transfer information. In addition to acoustic communication, visual communication also plays a major role in the animal kingdom. In an impressive way, this can be observed in marine-living cuttlefish (octopus). These visually communicating organisms are not only able to change their skin color and surface so that they are perfectly camouflaged. They are also able to move platelets in the chromatophores (specialized skin cells) to reflect specific patterns of polarized radiation. It is now thought that they have a language via these different skin changes, consisting of over 50 vocabularies, and can even form sentences in their specific sequence. Both acoustic and visual communication require that the sender and receiver have highly complex and specialized organs such as ears and eyes.

But what happens when, as in the case of bacteria, the organism cannot possess these organs because it consists of only a single cell? Although bacteria are only single-celled, they can communicate with conspecifics and bacteria of other species. This mechanism was first described in symbiotic bioluminescent bacteria by Nealson et. al. in 1970 and termed **Quorum Sensing** by Fuqua in 1994 and represents a biochemical language. The vocabulary of this language are so-called autoinducers, which are formed by the bacteria themselves and released into the surrounding space. The more bacteria exist in a spatial proximity to each other, the more of these signal molecules are present in this environment. When the concentration of the messenger molecules reaches a certain threshold value (the quorum), the autoinducers flow back into the bacteria, switch on some genes, switch off others and thereby abruptly change the behavior of the organisms. Sometimes the microorganisms form a biofilm that protects them from antibiotics and disinfectants. In other cases, they produce toxins and decomposing enzymes.

MSI® technology

Inspired by Nature | Addicted to Function

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Prof. Richard Hammond, Penrhos Bio Chief Executive Officer:
"We are thrilled to be working with our colleagues from **pro3dure medical** to combine the power of **Remora**™ technology with their world leading expertise in material development. The result is a completely novel way of inhibiting the formation of harmful multispecies biofilms in their products. **Remora**™ is inspired by the natural biology of certain seaweeds that are able to block communication between microorganisms. The result is a surface that remains clean and free of biofilms. Backed by Unilever, we have proven efficacy in a wide range of applications and we are really excited to bring this benefit into audio markets through the innovation of **pro3dure**, an ideal partnership."

MSI® – Award Winning Technology

Penrhos Bio | Unilever | pro3dure medical GmbH

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Winner of Scotland's Life Sciences Annual Awards 2024
Glasgow (Category Innovative Collaboration)
A great evening with our friends and partners was highlighted with a big win for our collaboration!

The formation of biofilms is therefore a significant medical problem. If the communication of microbes could now be interfered with by means of suitable molecules, bacterial infections could also be contained. Unlike antibiotics, the germs would not be killed, but cut off from their information channels. For this reason, the organisms cannot develop resistance to the interfering molecules, as is known from antibiotic resistance. **(Fig. 4)** In nature, biofilms are not observed on all moist surfaces. There are species of algae that are virtually free of biofilms. For example, the red alga Delisea produces a variety of natural compounds (furanones) that appear to block quorum sensing by bacteria. Inspired by nature, **Unilever / Penrhos Bio** has developed a biomimetic technology to combat biofilms (called **Remora**™). Research results also document that the bacteriostatic effect works not only on bacteria but also on yeasts, fungi and algae. Thus, a technology is available that can effectively combat „multi-species" biofilms.


In close cooperation between **Penrhos Bio** and **pro3dure medical**, this multi-species inhibition (**MSI® technology**) was realized for the first time in 3D printing and other audio materials. In doing so, **pro3dure medical** is taking an innovative step towards functional materials in the further development of its audio products in accordance with the **MSI®** guiding principle "Inspired by nature, addicted to function". Come and find out more details about **MSI® technology** and see numerous new, functional audio materials.

Abb. 3: **audioprint® 11.1 | MSI** mit **MSI® Technologie**.
Fig. 3: **audioprint® 11.1 | MSI** with **MSI® technology**.

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Since we have been using **audio**print® **GR-1** from pro3dure, we have achieved the best results in terms of break resistance, compatibility and precision.

Elmar Eck
Inhaber/owner LaborService Eck

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- exceptional UV stability
- high process robustness
- high elongation at break
- highest demands on aesthetics
- **MSI® technology** for the reduction of biofilms

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beige-opak / beige-opaque
Art.-Nr. / item-no.: **A1000860**

rötlich-transparent / reddish-transparent
Art.-Nr. / item-no.: **A1000863**

The generative manufacturing of earpieces on the base of digital data is nowadays a well established process in the HA business and represents an important part of different process chains. For this digital approach **pro3dure** medical offers a wide variety of technology packages and materials. Therefore, the light curing resins of the **audio**print® series cover the entire material spectrum of generative manufacturing. You can choose from resins for stereolithography to materials for image projection systems (DLP) with different irradiation sources. We offer customized material configurations on request. The **audio**print® **GR-1** is a light-curable polymerizable resin intended to be used in conjunction with curing light equipment, which is indicated for the fabrication of earmolds. The material distinguishes itself by excellent UV irradiation stability in combination with a high elongation at break. **audio**print® **GR-1** offers high process robustness and fulfills the highest demands on aesthetics without brittleness. The field of application for this high tech resin ranges from SLA technology to also DLP systems with UV sources ≤ 385 nm and irradiation sources ≤ 405 nm as **audio**print® **GR-10**.

In addition, we offer an **audio**print® **GR-1 | MSI** and **audio**print® **GR-10 | MSI** variant. The biomimetic 3D printing resins reduce, through the additional implementation of the **MSI® technology**, the adhesion of biofilms in a natural way to the surface of earmolds. This specifically improves hygiene in the ear.

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The new **audioprint® GR-11.1** is the 3D printing resin for the additive manufacturing of earmolds. In the further development of the proven **audioprint® GR-11**, the focus was on further optimizing the mechanical characteristics of the material in order to increase the stability of the 3D-printed earmolds and hearing aid shells once again. The exceptional mechanical stability compared to market competitors can be seen in **Figs. 5 to 7**. **audioprint® GR-11.1** earmolds have an almost 10-times higher total work of fracture, 5 times higher elongation at break, as well as a 3-times higher impact strength compared to the state of the art. In addition, we offer an **audioprint® GR-11.1 | MSI** variant. The biomimetic 3D printing resin **audioprint® GR-11.1 | MSI** reduces, through the additional implementation of **MSI® technology**, the adhesion of biofilms in a natural way to the surface of earmolds. This specifically improves hygiene in the ear.


Abb. 5: Pendelschlagfestigkeit (ASTM D256), J/m
Fig. 5: Pendulum impact resistance (ASTM D256), J/m

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- exceptional stability
- almost 10 times higher total work of fracture
- 5 times higher elongation at break
- almost 3 times higher impact strength
- extraordinary torsion resistance under cyclic fatigue loading (no damage even after 100 load cycles)
- TPO free material (no CMR[1] substances)
- **MSI® technology** for the reduction of biofilms

[1] carcinogenic, mutagenic, toxic for reproduction

audioprint® GR-11.1
1 kg Flasche / bottle
klar-transparent / clear-transparent
Art.-Nr. / item-no.: **A1001510**
schwarz / black
Art.-Nr. / item-no.: **A1001511**
hellbraun / light brown
Art.-Nr. / item-no.: **A1001512**
beige-opak / beige-opaque
Art.-Nr. / item-no.: **A1001513**
blau / blue
Art.-Nr. / item-no.: **A1001514**
rot / red
Art.-Nr. / item-no.: **A1001515**
weiß / white
Art.-Nr. / item-no.: **A1001516**

audioprint® GR-11.1 | MSI
1 kg Flasche / bottle
beige-opak / beige-opaque
Art.-Nr. / item-no.: **A1001520**
rötlich-transparent / reddish-transparent
Art.-Nr. / item-no.: **A1001523**

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The generative manufacturing resin **audioprint® GR-12.1** can be used particularly advantageous for the production of cast forms in order to realize a digital way for silicon earpieces. By defined diffraction particles a super smooth surface of the cast forms can be realized. The low viscosity, optimized mechanical properties and reduced inhibition layer helps you in the steps cleaning, post curing of the cast form and demolding the silicone earpiece.

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pro3dure´s audioprint® GR-12.2 generative resin has been specially developed for the casting process of polyurethane earmolds (e.g. **SysTherm®**). Due to the high green flex modus of the material thin cast molds can be realized. The slight translucency of **audioprint® GR-12.2** gives you the ability to observe the casting process in order to avoid trapped air in the final part. The material can be used in a favourable manner with the **SysTherm®** separator **SPT-1**.

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3D printing resins of the **audioprint® HPR series** are the ideal choice for additive manufacturing of hearing protection earmolds for use in work areas with the highest requirements. Hearing protection made of **audioprint® HPR-1** features outstanding chemical-physical properties such as the highest tensile strength in its class **(Fig. 8)**. In addition, **audioprint® HPR-1** earpieces have a reversible thermo effect. The material is "hard" when inserted into the ear canal and then becomes flexible due to body temperature, thus providing a comfortable and natural wearing comfort. This effect increases the acceptance of wearing hearing protection and thus occupational safety. Hearing protection made of **audioprint® HPR-1** can be optimally sealed with **audioprint® L-2 FLX** lacquer. **audioprint® HPR-1** is available in clear-transparent, white-opaque, blue-transparent, red-transparent and 3 neon signal colors.

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
Abb. 8: Reißfestigkeit, (ISO 37), MPa *
Fig. 8: Tensile strength (ISO 37), MPa *

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- outstanding properties
- highest tensile strength in its class
- reversible thermo effect
- natural wearing comfort
- increase of work safety

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audioprint® L-1 is the universal, radiation-curing one-component varnish for surface finishing of hard earmolds and custom shells (e.g. **audioprint® GR-1**, **GR-10** and **INR-1**). Its optimized flow and adhesion properties produce a homogeneous, high-gloss surface on the earpiece that is easy to clean and scratch-resistant. A novel additive system adapted to the acoustic application ensures maximum protection against yellowing, even after prolonged exposure to UV light. The varnish is available in a wide range of colors and is now also available as **audioprint® L-1 | MSI**, to reduce biofilm adhesion and facilitate cleaning of earmolds.

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- optimized flowing and adhesion characteristics
- **MSI® technology** for the reduction of biofilms

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The **audioprint® L-2 FLX** is a light-curing soft-touch coating for flexible and hard earmolds. Due to its novel chemical formulation, surfaces with special wearing comfort are achieved. The varnish adheres by both mechanical and chemical bonding, thus creating a long-lasting surface seal. **audioprint® L-2 FLX** is ideally suited for finishing generatively manufactured earmolds. Now also available as **audioprint® L-2 FLX | MSI** with reduced biofilm adhesion.

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- long-term surface sealant
- **MSI® technology** for the reduction of biofilms

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"Fischer Amps manufactures state-of-the-art in-ear headphones with the **audio**print® series from pro3dure. The quality of the materials is convincing and the support required for the complex processes is professional, uncomplicated and fast."

Jochen Fischer | Inhaber/owner Fischer Amps

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"DEER EAR IN EAR MONITORING's customers expect high quality and a wide range of designs. We can use the **audio**print® materials on a broad variety of printers in consistent quality and the wide range of colors helps us to implement our customers' design wishes and deliver high quality. High-quality sound meets high-quality materials."

Yogi Lang | Inhaber/owner DEER EAR IN EAR MONITORING

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- highest UV stability
- highest elongation values
- wide range of colors
- signal colors for industrial safety

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audioprint® **INR-1** generative resins have been developed for special demands like work protection and audio applications. All resins of the **audio**print® **INR-1 series** are characterized by highest UV-stability, extraordinary mechanical properties (e.g. highest elongation at break values) and a broad spectrum of colors. Four signal colors are available for 3D-printing of hearing protection earmolds and special colors for in-ear monitors like purple-opaque, blue-opaque and pink-opaque. Accordingly, the **audio**print® **INR-1** series is the perfect addition to the **audio**print® **GR-series** and offers you the possibility to fulfill your customer's demands more individual. **Make 3D-printing more colorful!**

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RPR-1 is the lightcuring thixotropic gel for the repair and modification of earpieces. The fine dispensing canulas guarantee optimal control and application of the material in order to make repair and venting works more convenient and efficient.

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- optimal control & application
- optimal post-processing

RPR-1	klar-transparent / clear-transparent	3 × 3,5 g	Art.-Nr. / item-no.: A1400100
RPR-1	rot-transparent / red-transparent	3 × 3,5 g	Art.-Nr. / item-no.: A1400200
RPR-1	blau-transparent / blue-transparent	3 × 3,5 g	Art.-Nr. / item-no.: A1400300
RPR-1	beige / beige	3 × 3,5 g	Art.-Nr. / item-no.: A1400400
RPR-1	Spritzenkanülen / syringes canulas	10 St. / 10 pc.	Art.-Nr. / item-no.: A4020134

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The light curing signature gel **SGNR-1** can be used to further individualize your earpieces. Due to its optimized thixotropy finest details of e.g. laser engraving can be filled without leaking or smearing during the working process. **SGNR-1** is available in red, blue and white. You can choose between a start packing (3 × 3.5 g in red, blue and white) or standard packing (2 × 3.5 g per color).

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- optimized thixotropy
- for finest details

SGNR-1	Startpackung / start package	3 × 3,5 g	Art.-Nr. / item-no.: A4020130
SGNR-1	Standardpackung rot / standard package red	2 × 3,5 g	Art.-Nr. / item-no.: A4020131
SGNR-1	Standardpackung blau / standard package blue	2 × 3,5 g	Art.-Nr. / item-no.: A4020132
SGNR-1	Standardpackung weiß / standard package white	2 × 3,5 g	Art.-Nr. / item-no.: A4020133
SGNR-1	Spritzenkanülen / syringes canulas	10 St. / 10 pc.	Art.-Nr. / item-no.: A4020134

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All under one roof

The 3D printing silicones for audiology applications from pro3dure medical are united under the umbrella brand otosil® to make your selection easier and clearer.

The possible applications for silicone materials in hearing aid acoustics are diverse and are constantly being expanded, for example, by combining them with 3D printing. This also increases the demands on the materials used. One example of this is the increasingly filigree design of silicone earmolds, which is only made possible by the casting process using 3D printing. This results in particularly high requirements with regard to e.g. the tear resistance of the earmold materials. Furthermore, the silicone varnishes used must have special flow and wetting properties in order to optimally reproduce the filigree details. All materials in the **otosil® series** have been designed and developed with these parameters in mind. **Now with MSI® technology!**

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otosil® IS-2 with integrated color indicator

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otosil® IS-2 is our new universal impression material with integrated color indicator from ocean-blue to white. In this way, the setting process of the material in the patient's ear can be perfectly monitored by the acoustician. Furthermore, **otosil® IS-2** is characterized by highest quality, creamy initial consistency as well as optimized thixotropy. The squeeze properties have been optimized even further. The material is sharp to the detail, has an adjusted tear resistance, a good recovery from deformation, is scannable, easy to process, skin-friendly and biocompatible (A-silicone, medium viscosity, ear retention time approx. 2-3 min., final hardness approx. 35 Shore A).

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- integrated color indicator
- controlled setting process
- easy to squeeze out due to optimized squeeze properties
- well scannable
- biocompatible



Abb. 9: Dosierkraft (Spiralmischer S50), N
Fig. 9: Dispensing force (spiral mixer S50), N

* Gemäß den internen Designvorgaben.
According to internal design specifications.

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otosil® PS-1 is the universal impression material with a soft-creamy consistency. The soft, non-sticky putty material is particularly easy to mix. Optimum mixing is easily recognized by the orange/white color combination. The material is characterized by perfect sharpness of detail, has excellent tear resistance and sets quickly. **otosil® PS-1** is scannable, easy to process, skin-friendly and biocompatible (A-silicone, kneadable consistency, ear setting time approx. 3-5 min., final hardness approx. 40 Shore A).

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- soft-creamy consistency
- non-sticky
- easy to mix
- high detail sharpness
- scannable
- biocompatible

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- optimized hardness (50 Shore A)
- 3 signal colors
- excellent damping properties
- **MSI® technology** for the reduction of biofilms

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otosil® protect 50 | MSI is our material for durable silicone-based hearing protection. The material offers various advantages such as excellent adaptability and high wearing comfort. Thanks to the specifically optimized hardness of 50 Shore A, both very soft and more stable earmolds can be manufactured with just one material. The excellent mechanical properties enable both particularly thin-walled earmolds and conventional geometries. At the same time, the earmolds are characterized by significant robustness and resistance. With **MSI® technology** for improved hygiene.

otosil® protect 50 | MSI
8 x 48 ml Doppelkartuschen à 2 x 24 ml A+B /
8 x 48 ml double cartridges of 2 x 24 ml A+B

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The earmold silicones of the **otosil® ES-1.1 series** offers everything you expect from a modern high-strength silicone material for the production of flexible earpieces. Available in Shore A hardnesses of 25, 40 and 60 and a broad spectrum of different colors. Due to the optimized flow characteristics and low dispensing pressure generative fabricated cast forms can be filled without air traps in an efficient and controllable way. By a newly developed crosslinker matrix an extraordinary combination of high elongation and high tear and tensile strength is achieved. This enables you to remove the earpiece from the cast with a reduced tendency of microcuts. The working time in a pressure pot is at 50 °C, 60 min. Now also available as **otosil® ES-1.1 | MSI** with reduced biofilm adhesion.

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- broad spectrum of colors
- optimized flow characteristics
- high efficiency and control
- combination of high elongation, high tear and tensile strength
- **MSI® technology** for the reduction of biofilms

8 x 48 ml Doppelkartuschen à 2 x 24 ml A+B / 8 x 48 ml double cartridges of 2 x 24 ml A+B

Farbe / color		25 Shore A Art.-Nr. / item-no.	40 Shore A Art.-Nr. / item-no.	60 Shore A Art.-Nr. / item-no.
otosil® ES-1.1	klar-transparent / clear-transparent	A0030800	A0032200	A0033600
	rötlich-transparent / reddish-transparent	A0030203	A0031603	A0033003
	schwarz-opak / black-opaque	A0030376	A0031776	A0033176
	weiß-opak / white-opaque	A0030100	A0031500	A0032900
	neon-grün-opak / neon-green-opaque	A0030530	A0031930	A0033330
	neon-gelb-opak / neon-yellow-opaque	A0030729	A0032129	A0033529
	neon-pink-opak / neon-pink-opaque	A0030903	A0032303	A0033703
	neon-orange-opak / neon-orange-opaque	A0031327	A0032727	A0034127
	dunkelblau-opak / dark blue-opaque	A0030428	A0031828	A0033228
	hellblau-opak / light blue-opaque	A0031128	A0032528	A0033928
	rot-opak / red-opaque	A0030627	A0032027	A0033427
	violett / violet	A0030680	A0031680	A0033004
	rauchschwarz / smoke grey	A0030690	A0031690	A0033690
otosil® ES-1.1 \| MSI	neon-gelb-opak / neon-yellow-opaque	A0030750	A0032130	A0033010
	neon-pink-opak / neon-pink-opaque	A0030751	A0032131	A0033011
	neon-orange-opak / neon-orange-opaque	A0030752	A0032132	A0033012

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one-component lacquer for sealing silicone earmolds

otosil® SL-1 is ideal for smoothing irregularities, such as those that can be transferred to the surface in the "digital workflow" by building steps in 3D-printed molds. The varnish is characterized by optimized flow and wetting properties in combination with exceptional adhesion to silicone surfaces. In addition, it is free of aromatic solvents and thus supports you in the laboratory in terms of occupational safety. Otoplastics coated with **otosil® SL-1** literally „shine" with their surface quality. Now also available as **otosil® SL-1 | MSI**, to reduce biofilm adhesion and thus facilitate the cleaning and hygiene of earmolds.

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- optimal smoothing and easy processing
- optimized flow and wetting properties
- exceptional adhesion to silicone
- fast setting (30 min. at 80 °C)
- no aromatic solvents
- no unpleasant odor
- **MSI® technology** for the reduction of biofilms

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

For the production of polyurethane-based thermoelastic earmolds, **SysTherm®** from pro3dure has been our top priority since 2019. **SysTherm®** has made our production faster and more reliable.

Salvatore Falcone
Inhaber/owner Otoprint GmbH

Abb. 11 (a-j): Herstellungsschritte **SysTherm®**. / Fig. 11 (a-j): Production steps **SysTherm®**.

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

- no toxic labelling
- available in 2:1 cartridges
- easy to use & less waste
- linear thermo effect
- highest elongation at break
- superior mechanical properties
- **training with certificate**

SysTherm® PU
Standard-Packung / standard package
2 x 50 ml + 4 Mischkanülen /
2 x 50 ml + 4 mixing canulas

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thermo effect – no toxic ingredients

From now on you work systhermatically! According to the holistic approach of the company **pro3dure**, we offer you with the product family **SysTherm®** everything to produce earmolds with thermo effect efficiently and under highest quality standards. "From now on you work systhermatically" is to be taken literally.

The **SysTherm®** system consists of a unique 2-component material for earmolds with thermo effect. **SysTherm®** combines two world firsts. In addition to outstanding chemical and physical properties, the **SysTherm®** plastic has no toxic labeling and is also available in a 2:1 cartridge version. This means that not only the end customer receives an earmold with outstanding wearing comfort, but also the laboratory employee is assured of maximum production efficiency and occupational safety. **SysTherm®** is rounded off by the **SysTherm® SPT-1** separator, which has been optimized for **SysTherm®** plastic and can be used with conventional 3D printing resins for cast molds (e.g. **audioprint® GR 12.2**). This multitude of innovations of the **SysTherm®** system thus increase the efficiency and robustness of your laboratory processes with maximum flexibility.

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PR-1 is a light-curing pouring resin for the traditional production of earmolds, IE and IEM shells. High elongation at break and excellent UV-stability in combination with **pro3dure´s** curing device **CD-2** comply with all demands on modern light curing technology. This provides user-friendliness and cost effectiveness. The material is available in a broad variety of colors.

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The water-soluble **CL-1** cleaning agent removes resin residue from additive molded parts more effectively and gently than conventional Isopropanol.

- 5 times higher cleaning effect°
- 5 times higher cleaning capacity°
- non-explosive / non-flammable
- ideal for ultrasonic systems
- highly efficient
- ultrafiltrable
- compatible with cleaning devices such as the **CLD-1**
- low-odor

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For your individual needs different scanner models are available:

- **free Stage one:** Stylish and impressive - the elegant solution for moderate production requirements in the store.
- **free Stage twin:** Efficient, versatile and future-oriented - the complete package for store use.
- **mono Scan:** The reliable entry-level solution for hearing care professionals - high precision at a low price.
- **duo Scan:** Proven technology that pays off - the classic for fast, high-quality scans.
- **ds Production:** Robust, compact and precise - the marathon runner for continuous use under production conditions.

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
An important building block for the digital workflow is the ear impression scanner. We work in close collaboration with our partner smart optics in order to achieve flexible solutions for your laboratory or shop. Various scanner models are available for this purpose, all of which are based on the proven stripe light triangulation.

The **free Stage twin** was developed for quality-conscious audiologists and is the ideal solution for trend-oriented shops. With its slim design, the system-independent 3D scanner is easy to set up at any workstation. The smooth surfaces as well as the wide opening allow easy access. The scanner's ergonomic shape, obtained by deliberately not using a cover for the housing, simplifies the working processes and the ultra-quiet motors give the user the opportunity to let his customer experience the entire scanning process impressively. The innovative touch sensor makes it possible to start the scanning process directly on the scanner.

This saves you valuable time. The consistently high quality of materials used in smart optics products is always most important.

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Kralle / Claw

CAD/CAM Software für Profis

CAD/CAM software for professionals

Fully in line with our process orientated company philosophy and in order to complete the digital workflow we offer you the CAD/CAM packages of our partner **Cyfex**.

Good earmolds provide a service inconspicuously. Through their shape and their material, they fulfil multiple requirements - acoustics and comfort, among others. You can rely on the **Secret Ear System** with regards to quality and the satisfaction of your clients. **Secret Ear Designer** is the high-performance earmold modeling software for the digital manufacture of ear protection, earmolds and in-ear-monitoring systems. The application facilitates a comfortable, intuitive and time-saving analysis for users and distinguishes itself through a high degree of automation. With only one license, BTE as well as ITE earmolds can be modeled. Digital technology offers you higher quality and thus less and uncomplicated remakes. Thanks to high performance, you can obtain higher output with the same number of employees. In addition, you have the advantage of being able to develop your own products and have a tool at your disposal to create your own construction data.

- **All earmold styles**
 With just one Secret Ear Designer license you can design IdO and HdO earmolds, as well as ear protection and in-ear-monitors, including molds for soft earmolds.
- **Intuitive interface**
 Thanks to easier handling, employees are in a position to model earmolds that they have designed independently and in detail after a few hours of training.
- **Time savings**
 On the basis of easily adjustable product templates the required modeling steps can be reduced to an absolute minimum.
- **ERP connectivity**
 Two-way transfer of parameters between order management system and Secret Ear Designer by means of ERP-Link.
- **Compatibility with output media**
 The application is compatible with all 3D printing technologies available on the market as well as with milling centers for subtractive production.
- **Flexible shape design**
 For maximum independence and creative freedom, construction elements can be autonomously imported and/or designed.

Secret Ear Builder provides your lab with a tool for complete print job generation. The intuitive software makes maximum use of the printer's build area and allows fast grouping and automated alignment of the earmolds to be printed. All necessary algorithms are available for the creation of support structures optimized for earmolds. All relevant basic settings – also those of your 3D printer – are preconfigured.

Different attractive license models are available for this **CAD/CAM package**. Let us find the best fitting solution for your individual demands.


Download-Link SysTherm® Cyfex-Vorlagen
https://www.pro3dure.com/de/cat/index/sCategory/121

Download link SysTherm® Cyfex templates
https://www.pro3dure.com/en/cat/index/sCategory/121

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The new **ASIGA MAX 2™** forms the backbone of your digital process chain with the **audioprint®** materials. This 3D printer, specially developed for the digital workflow in the laboratory, is characterized by maximum productivity with minimum space requirements. With a pixel size (x,y) of 62 µm and a 385 nm UV irradiation source, you have access to the complete material portfolio of the **pro3dure audioprint® series**. This makes the ASIGA MAX 2™ 3D printer an all-rounder.

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Discover the new **Asiga Ultra™**. The Active Layer Monitoring technology guarantees precise layer thickness control for maximum print accuracy and speed. Experience razor-sharp 4K image resolution over a large build area. With the auto-calibrating UV-LED system and DLP micromirror technology, the Asiga Ultra offers reliability and versatility. Touchless entry allows you to interact effortlessly without contaminating surfaces. The magnetic platform holder and infrared material heater optimize your workflow and enable the use of a wide range of advanced **pro3dure audioprint®** materials. The intuitive user interface and Transparent Mode™ make the printing process easier than ever. Rely on chamber heating for optimal process reproducibility.

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The new **MiiCraft Alpha 2 K** desktop 3D printer is characterized by its compact size and offers ideal application possibilities for audiology. With a build space of up to 130 mm in length and exceptionally high resolution for precise printing (1980 x 1200 pixels with XPR™), it enables extremely accurate results. The build platform and resin tank are designed for user-friendly handling and the acrylic covers are open on three sides to provide a larger workspace.

RFID (material identification)
- traceability (resin batches, print platform, resin consumption)

High durability
- LED life reaches up to 20000 hours
- tank life up to > 7.5 kg

- excellent sharpness throughout the image
- sharpness MTF > 50 %, distortion < 0.1 %

High speed – high precision
- platform and tank easily removable
- 3-sided open acrylic cover for spacious work area
- increased print and file transfer speed

Technische Daten | technical data

Technische Daten / technical data	Miicraft Alpha 2K
Bau-Volumen (X, Y, Z) / build volume (X, Y, Z)	130 x 73 x 75 mm
Laterale Auflösung / lateral pixel resolution	68 µm
Z-Layer Auflösung / Z-layer resolution	5-500 µm
Technologie / technology	DLP
LED-Wellenlänge / LED wavelength	385 nm
User interface	5 Zoll Touch Screen / 5 inch touch screen
Dimensionen / system sizing	30 x 40 x 42,26 cm / 16 kg
Input/Output	Ethernet, USB, WIFI Dongle

Miicraft Alpha 2K
Art.-Nr. / item-no.: **U2017005**

CLD-1 Cleaning Device

CLD-1
Cleaning device / cleaning device
Art.-Nr. / item-no.: **U2020011**

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The "cleaning device" **CLD-1** represents another building block for the optimized and controlled manufacturing process of medical devices by generative fabrication. Unlike conventional ultrasonic baths the **CLD-1** generates a vortex of the cleaning agent. As a result, microcracks in the component surfaces are significantly reduced, thus further increasing the part quality in terms of surface finish, appearance and mechanics of your objects. A selection of different cleaning programs assists you in perfectly matching the cleaning program to the selected material class and choice of cleaning agent. The cleaning container can be easily removed and changed in order to simplify the manufacturing processes in your lab. Changing the cleaning liquid can be done within seconds in a convenient way.

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The **RLR-1** is a roller mixer for homogenization of 3D printing materials with pigment or filler content. 3D printing materials like the **audioprint® INR-1** require a thorough mixing before use, handshaking is often slow. This can result in significant colour deviations and even build failures may occur. The **RLR-1** Mixer is constructed for constant operation in order to keep your **pro3dure** 3D materials ready for use at any time with the same consistency.

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
The high tech post-curing unit **CD-2** is tuned to the latest generation of 3D printers. It features top-level productivity, the most current LED irradiation technology, sophisticated sensors and a wide range of smart features. Through the touch panel you can fully control and monitor all relevant irradiation parameters. Sensors measure the protective gas flow and inform the user about the status of the light source. This avoids inhibition layers on generated parts. As a result, the biocompatibility and the optimal mechanical properties of your 3D printed parts are guaranteed. In a profound and unprecedented form the **CD-2** combines compactness and efficiency on the highest level of technology, along with an excellent price/performance ratio.

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USA
pro3dure medical LLC
Valley View Business Center
9825 Valley View Road
Eden Prairie, MN 55344
Phone: +1 (952) 426 1928

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Exclusion of liability: We cannot be held responsible for false statements in this brochure.
Delivery of our products is effected according to our general terms of delivery and payment. Technical modifications reserved.